



Asahi**KASEI**
BIOPROCESS

Understanding Planova™ Virus Filtration and Designing Robust Virus Filtration Processes


Tomoko Hongo-Hirasaki, Ph.D.
Senior Scientist
Group Manager of Evaluation and Analysis
Asahi Kasei Medical Co., Ltd.

Nov 2017, Prague, 20th Planova Workshop

1. Introduction
2. Findings for Planova Virus Filtration
 - 2.1. Key Factors of Robust Virus Filtration Processes
 - 2.2. Design Space for Robust and Stable Performance
3. Findings for Appropriate Viral Clearance Studies
4. Conclusions

Line-Up of Planova Products

- ❑ Pioneering nanofiltration technology since 1989



Planova Filters



Filter	Planova N	Planova BioEX
Hollow fibers	Regenerated cellulose	PVDF
Pore size (nm)	15, 20, 35, (75)	Parvovirus removal
Filtration surface (m ²)	4.0, 1.0, 0.3, 0.12, 0.01, 0.001	4.0, 1.0, 0.1, 0.01, 0.001, 0.0003
Max TMP (kPa)	98	343
Water flux (LMH)	20N: 66 @ 98 kPa	170 @ 343 kPa
SIP (Autoclaving)	NO	YES
Post-use integrity test	Leakage test + Gold particle test	Leakage test

- Viruses bigger than the pore size of the filter are removed.





Size exclusion mechanism

- Voids are connected 3-dimensionally with several capillaries.
- This structure is repeated through the membrane thickness.




Multi-layer filtration mechanism




Negative image

Multi-layer image




Mechanism: Size Exclusion

Correlation between virus size and removal on Planova 35N



Designed Pore Size Gradient Membrane









Planova 20N membrane
Model gold particles: 10, 20, 30 nm
Microscopic analysis

- 30 nm gold particle trapped at inner side as a function of pre-filter layer
- 20 nm gold particle trapped at center of membrane

⇒ High capacity


- Resistance to plugging
- Stable flux performance
- High virus trapping capability

Parvovirus Trapping Images




5 mg/mL IgG/100 mM NaCl (pH 4.5)

Flux scalability



Viral clearance scalability



- Robust performance and viral clearance was observed for Planova 20N filters at all scales

Applicable to Wide Range of Proteins

Pharmaceutical	MW (kDa)	P15N	P20N/PBioEX	P35N
Interleukin	12 – 35			
Interferon, GCS-F	20			
MCS-F, EPO, Thrombin	30 – 40			
Urokinase	54			
α 1-protease inhibitor	56			
Factor-IX	58			
Protein C	62			
Hemoglobin, antithrombin	65			
Albumin	66			
Osteopontin	70			
Monoclonal IgG, Polyclonal IgG	160			
Factor XI	180			
Factor VIII	300			
Fibrinogen	340			
IgM	800			

Involvement of Virus Filter in Drug Development

Asahi**KASEI**
BIOPROCESS



Easy optimization

- Stable filtration
- Reliable virus removal
- Time saving
- Multi-product applicability

Easy to platform

Scalability

- Reliable scalability
- Reliable virus removal
- Time and cost saving

Consistency

- Stable and reliable filtration
- Assured virus safety
- Time and cost saving

Robust performance

Most important feature of virus filter

1. Introduction
- 2. Findings for Planova Virus Filtration**
 - 2.1. Key Factors of Robust Virus Filtration Processes**
 - 2.2. Design Space for Robust and Stable Performance**
3. Findings for Appropriate Viral Clearance Studies
4. Conclusions

1. Introduction
2. Findings for Planova Virus Filtration
 - 2.1. Key Factors of Robust Virus Filtration Processes**
 - 2.2. Design Space for Robust and Stable Performance
3. Findings for Appropriate Viral Clearance Studies
4. Conclusions

- ✓ Filterability
 - Stable flux
 - Small flux decay
 - Predictable throughput filtration to target processing time
 - Scalability; consistent performance across filter sizes



- ✓ Virus Removal
 - Maintain virus LRV >4
 - Reliable virus removal rate
 - Scalability; consistent performance across filter sizes

Factors Necessary for Robust Virus Filtration

- ✓ Process parameters
 - Protein conc.
 - Solution property
pH, ionic strength etc.
 - Temperature
 - Filtration pressure (flow rate)

- ✓ Filtration scale
 - Lab to manufacturing
Scalability, consistency

What is Robust Filterability?




Pressure: 0.8 bar (P20N), 3.0 bar (PBioEX)

Product solution: 10 or 30 mg/mL IgG, 100 mM NaCl, pH 4.5

- ✓ Planova filters show stable flux with relatively small flux decay and predictable throughput for the duration of the filtration process.

What is Robust Filterability: Impact of IgG Concentration on Flux Decay


Product solution: 10 or 30 mg/mL IgG, 100 mM NaCl, pH 4 – 4.6




- ✓ Planova BioEX shows low flux decay even for mid to high IgG concentration.

What is Robust Filterability: Impact of pH and ionic strength on flux decay


Planova 20N



Planova BioEX



Competing Filter A



10 mg/mL human IgG, 3 h filtration

Relative flux (%)
$$= (3\text{h flux} / \text{Initial flux}) \times 100$$


- ✓ Planova 20N and BioEX show stable flux over wide range of pH and ionic strength

What is Robust Filterability: Impact of Multimer on Flux Decay

Pressure: 3.0 bar

Product solution: 10 mg/mL IgG, 100 mM NaCl, pH 4.5,
0.0025% multimer spiked (low pH denaturation)

Prefilter: Planova™35N (post-multimer spiking)



- ✓ Planova BioEX is less impacted by aggregates and shows lower flux decay compared to Filter A.
- ✓ Multimer caused by changes in parameters may impact filterability

Robust Virus Filtration: Impact of Protein Concentration on LRV

AsahiKASEI
BIOPROCESS

 **Parvovirus LRV
(Antibody Concentration)-20N**

IgG Conc. (mg/ml)	PPV LRV	Filtration Time (hr)	Filtration Volume (L/m ²)	Throughput (kg/m ²)
1	≥ 5.67	6.4	350	0.35
5	≥ 5.37	6	280	1.40
10	≥ 6.00	6	240	2.40
20	≥ 5.50	6	148	2.96
30	≥ 5.58	6	98	2.94
50	≥ 5.67	5	24	1.20

100 mM NaCl, pH 4.5,
Pressure; 78.4 kPa,
Serum-free PPV 0.5 vol%

Connecting People, Science and Regulation

 **Parvovirus LRV
(Antibody Concentration)-BioEX**

IgG Conc. (mg/ml)	PPV LRV	Filtration Time (hr)	Filtration Volume (L/m ²)	Throughput (kg/m ²)
1	≥ 5.42	4	490	0.49
5	≥ 5.78	4	470	2.35
10	≥ 5.35	4	425	4.25
30	≥ 5.28	6.4	300	9.00
50	≥ 5.10	9.4	180	9.00

100 mM NaCl, pH 4.5
Pressure; 294 kPa,
serum-free PPV 0.5 vol%

Connecting People, Science and Regulation

Tomoko Hongo-Hirasaki, 2011 PDA Virus & TSE Safety Conference (adapted)

- ✓ No impact on PPV LRV (>4 log) over wide range of protein concentrations

Robust Virus Filtration: Impact of pH on LRV

 **Parvovirus LRV (pH) - 20N**

pH	PPV LRV	Filtration Time (hr)	Filtration Vol. (L/m ²)
4.5	≥ 6.00	6.0	240
5	≥ 6.00	6.0	232
6	≥ 6.00	6.0	217
7	≥ 5.78	6.0	205
8	≥ 6.09	6.0	213

10 mg/ml IgG, 100 mM NaCl
Serum-free PPV 0.5 vol%,
Pressure; 78.4 kPa

Connecting People, Science and Regulation®

 **Parvovirus LRV (pH) - BioEX**

pH	PPV LRV	Filtration Time (hr)	Filtration Vol. (L/m ²)
4.5	≥ 5.28	6.4	300
5.6	≥ 5.56	9.3	215
6.7	≥ 5.14	9.0	150
7.6	≥ 5.40	9.6	155

30 mg/ml IgG, 100 mM NaCl
Serum-free PPV 0.5 vol%,
Pressure; 294 kPa

Connecting People, Science and Regulation®

Tomoko Hongo-Hirasaki, 2011 PDA Virus & TSE Safety Conference (adapted)

- ✓ No impact on PPV LRV (>4 log) over wide range of pH

Robust Virus Filtration: Impact of Ionic Strength on LRV

 Parvovirus LRV (Ionic Strength) - 20N

NaCl Conc. (mM)	PPV LRV	Filtration Time (hr)	Filtration Volume (L/m ²)
1	≥ 5.84	6.0	152
10	≥ 6.34	6.0	200
100	≥ 6.00	6.0	240
200	≥ 5.67	6.0	220
500	≥ 6.00	6.0	190

10 mg/ml IgG, pH 4.5
Serum-free PPV 0.5 vol%,
Pressure; 78.4 kPa

Connecting People, Science and Regulation® 19

 Parvovirus LRV (Ionic Strength) - BioEX

NaCl Conc. (mM)	PPV LRV	Filtration Time (hr)	Filtration Vol. (L/m ²)
50	≥ 5.48	6.0	300
100	≥ 5.28	6.0	300
250	≥ 5.28	9.5	280
500	≥ 5.92	9.5	160



30 mg/ml IgG, pH 4.5
Serum-free PPV 0.5 vol%,
Pressure; 294 kPa

Connecting People, Science and Regulation® 20

Tomoko Hongo-Hirasaki, 2011 PDA Virus & TSE Safety Conference (adapted)

- ✓ No impact on PPV LRV (>4 log) over wide range of ionic strength

Robust Virus Filtration: Impact of Temperature on LRV




10 mg/ml IgG, 100 mM NaCl, pH 4.6
Serum-free PPV, 0.5 vol%
Pressure: 98 kPa

20 mg/ml IgG, 100 mM NaCl, pH 4.5
Serum-free PPV, 0.5 vol%
Pressure: 343 kPa

- ✓ No impact on PPV LRV (>4 log) over wide range of temperature

Impact of Filtration Pressure

- ❑ The effects of low pressure and depressurization are dependent on other factors



Other factors:

- pH
- NaCl concentration
- Product
- Virus challenge

2011 PDA Virus and TSE Safety, Tomoko Hongo-Hirasaki, Asahi Kasei Medical

Process Pause Impact


- ❑ In 2011, the impact of pressure release was presented at the PDA Virus and TSE Safety Meeting

Virus Breakthrough after Pressure Release during Virus Retentive Filtration

27-Jun-2011 PDA, Barcelona
Dr. Marcel Asper

charles river

► pressure release



charles river

Most filter types display significant impact on parvovirus throughput due to the process pause

- Based on size exclusion mechanism, Planova filters demonstrated robust performance across broad factors;
protein conc., ionic strength, pH, temperature

- Pressure control is very important for virus filtration.






We investigated:

- ✓ What is the additional mechanism to explain these phenomena?
- ✓ Design space study for robust and stable performance across broad ranges of parameters

1. Introduction
2. Findings for Planova Virus Filtration
 - 2.1. Key Factors of Robust Virus Filtration Processes
 - 2.2. Design Space for Robust and Stable Performance**
3. Findings for Appropriate Viral Clearance Studies
4. Conclusions

Visualization of Virus Distribution Inside Membrane

Cross-sectional view of membrane using fluorescence microscopy




- ✓ Small virus entrapment location may depend on factors such as solution conditions and filtration pressure.

- ❑ Physicochemical properties (solution condition and/or molecule type) may impact the size-exclusion mechanism

Two hypothetical cases

Pore size reduction




Pore size reduction by adsorption of protein on membrane

Increase in virus removal capability

Decrease in flux

Aggregation



Increase in interactions, complexes act as large particles

Increase in virus removal capability

Decrease in flux


- ✓ Possible interactions: surface charge, or hydrophobicity of virus, protein and membrane

Adapted from Tomoko Hongo-Hirasaki, 2013 PDA Virus & TSE safety

Additional Mechanism: Mechanistic Effect

❑ Constraint by hydrodynamic force


Hydrodynamic force controls Brownian motion of viruses



u , Flow velocity is controlled by filtration pressure
 U_{cr} Critical velocity to overcome Brownian motion

Mechanism 1: Size exclusion $d_{virus} > d_{cap}$
 Mechanism 2: Hydrodynamic force

❑ Performance and process parameters





- Study goals: Define the robust operating spaces for Planova filters
 - ✓ Design spaces for Planova 20N and BioEX filters which provide highly effective parvovirus removal ($LRV \geq 4$)


Design Space Study: Design Space of PPV LRV

Human IgG

20N



BioEX




- ✓ Planova 20N and BioEX showed robust performance over wide range of conditions in design space studies.




Collaboration with Janssen

- Three mAb products (IgG1)
- Spiked with MVM ($\sim 6 \log \text{TCID}_{50}/\text{mL}$)

Filtration scheme for initial scouting and design space




Filtration scheme for low pressure limit determination



For Planova 20N and BioEX, operating pressures greater than 7.3 and 10.0 psi (0.50 and 0.69 bar), respectively, are expected to provide LRV greater than 4 for most mAb processes

Implementation in DSP





- ✓ Flexible process optimization possible with placement of virus filtration step
- ✓ Platform possible with applicability to wide range of conditions

1. Introduction
2. Findings for Planova Virus Filtration
 - 2.1. Key Factors of Robust Virus Filtration Processes
 - 2.2. Design Space for Robust and Stable Performance
- 3. Findings for Appropriate Viral Clearance Studies**
4. Conclusions

Considerations for Viral Clearance Studies

❑ Procedure for viral clearance study





Impacts of Overspiking Virus

❑ CRO experience:

Protein /run	Logs virus loaded (PFU)	LRV filtrate pool	Comments
A1	9.2	≥7.6	Inconsistent clearance between replicates
A2	9.3	6.1	
B1	9.2	≥7.4	
B2	9.1	5.8	
C1	9.4	3.2	Low overall reduction
C2	9.5	3.3	
D1	9.0	3.1	
D2	8.9	3.4	

Observation



- On filters from multiple suppliers, overspiking can lead to inconsistent results between duplicate runs or low parvovirus clearance

Viral Clearance Studies: Impact of Virus Load on LRV

Planova 20N

10 mg/mL h-IgG; 0.8 bar; pH 7, 100 mM NaCl



Virus spike: 0.03, 0.3, 3.0 vol%



✓ Pre-pause PPV LRV increases with virus load.

Viral Clearance Studies: Impact of Virus Load on LRV

Planova 20N
10 mg/mL h-IgG; 0.8 bar; pH 7, 100 mM NaCl
Virus spike: 0.03, 0.3, 3.0 vol%



- ✓ Pre-pause PPV LRV increases with virus load.
- ✓ It is important to consider the optimal balance of target LRV and virus load.
- ✓ Select appropriate virus load to obtain the target LRV.

- The virus load solution contains a high number of non-infectious virus particle content in addition to infectious particles.
- The total virus particle number could be more than expected and it would sometimes cause unreliable LRV.
- Ratio of non-infectious to infectious particles is variable for parvovirus (e.g., MVM): c.a., 250:1 – 2,500:1

Sampling of references:

- ✓ *Virology*, 2015, 481: 63-72
 - ✓ *J. Virol.*, 2005, 79: 289-98
 - ✓ *J. Virol.*, 2010, 84(15):7782-92
- In consideration of the limited virus capture capacity of virus filters, selection of virus load sufficient to achieve >4 log reduction is important.

❑ Select appropriate virus load in virus filtration studies


- Spiking appropriately can minimize assay artifacts
 - We recommend not to use spiking based on percentage ratio
 - Spike to a targeted total viral challenge
 - Typically, spiking to a total virus challenge of 10.5-11.0 log IU/m² is sufficient to obtain a virus LRV of 5-6
- 7.5–8.0 log IU on a 10 cm² filter

Example:

Method	1% spike ratio	Target challenge 8.0 log IU
Load volume	300 mL	300 mL
Virus stock titer	9 log IU/mL	9 log IU/mL
Spike volume	3 mL	0.1 mL
Viral load	7.0 log IU /mL	5.5 log IU/mL
Total viral load / 10cm ²	9.5 log IU	8.0 log IU

Considerations for Viral Clearance Studies

□ Procedure for viral clearance study



Points

- ✓ Select appropriate virus load to obtain the target LRV

- ✓ Plan at least the same process pause duration as in manufacturing process (e.g., for post-wash)

* Titration

1. Introduction
2. Findings for Planova Virus Filtration
 - 2.1. Key Factors of Robust Virus Filtration Processes
 - 2.2. Design Space for Robust and Stable Performance
3. Findings for Appropriate Viral Clearance Studies
- 4. Conclusions**

❑ To develop robust virus filtration processes

1. Planova virus filtration

- ✓ Based on size exclusion mechanism, Planova filters demonstrate robust small virus removal and filterability across broad conditions.
- ✓ Other mechanisms (factors) must be considered on a case-by-case basis (each product and process) in order to achieve higher small virus removal capability.
- ✓ Planova filters shows robust performance over wide range of conditions in design space studies.

2. Evaluation by appropriately designed viral clearance study

Selection of appropriate virus load and experimental design can improve outcomes and avoid study artifacts in viral clearance studies.

- **Asahi Kasei Medical Co., Ltd.**

Marcus Inouye
Yoshiro Yokoyama
Kouichiro Yanagida
Takuma Iwasaki
Ryota Oguri
Hirohisa Hayashida
Masayasu Komuro
Ryo Hamamoto
Akika Futamura
Yoshiyuki Sawamura
Yuma Ueno
Akiyo Kajitani
Sachiko Goshi
Mana Nakashima

- **Janssen**

Josh Goldstein
Dominick Vacante

- **WuXi AppTec, Inc.**

Michael Burnham
Joe Hughes

- **Asahi Kasei Bioprocess America, Inc.**

Daniel Strauss
Kazuya Kobayashi
Brian Buesing

- **Asahi Kasei Co., Ltd.**

Aguru Yamamoto
Yoshihiko Uchi