



PennState




Flow Distribution in Commercial Depth Filter Capsules Performance and Scale-up

Negin Nejati, Minyoung Kim, Andrew Zydney
Lara Fernandez Cerezo, Ehsan Espah, David Roush

November 8th, 2021


Depth Filtration in Monoclonal Antibody Processing

Depth filters used for both clarification and prefiltration



The global depth filtration market is projected to reach \$2.9 Billion by 2025

Depth Filter Operations



- Species removed by a variety of mechanisms including
 - Surface retention
 - Straining (within filter)
 - Adsorption
 - Settling
- Currently: Empirical sizing with safety factors for Bulk Drug Substance

Roush DJ, Lu Y. 2008. Advances in Primary Recovery: Centrifugation and Membrane Technology. *Biotechnol. Prog.* **24**:488–495.
<http://doi.wiley.com/10.1021/bp070414x>.

Depth Filter Scale-Up

Clarification of 2000L bioreactor using 0.25 m² depth filter

- Using the manufacture recommendation and safety factor of 1.5
- **225% excess footprint!**

Supracap™ 50
Small-scale



Scale-up


Stax™
Large-scale



Surface area per capsule:
 $0.22 \times 10^{-4} \text{ m}^2$

Surface area per capsule:
 0.25 m^2

TMP Profiles



Objectives

Develop a fundamental understanding of factors governing performance and scale-up of depth filtration processes using Pall PDH4 lenticular media as model system




1. Develop CFD model describing flow profiles in scale-down and pilot-scale depth filter modules
2. Validate CFD model using dye-binding and residence time distribution measurements
3. Evaluate DNA removal characteristics in both scale-down and pilot-scale modules

Model Depth Filter – PDH4 Supracap™ 50





PDH4 Supracap™ 50
4-Layers of depth
filter

Cellulose Fibers





2 mm

Diatomaceous
Earth (DE)




SEM image of the media cross section

A binder



CFD Analysis

PDH4 capsule
22 cm² surface area



Depth Filter Media

$$\rho(\mathbf{u}^* \cdot \nabla \mathbf{u}^*) = -\nabla p + \mu \nabla^2 \mathbf{u}^* - \frac{\mu}{\kappa} \mathbf{u}^*$$

$$\nabla \cdot \mathbf{u}^* = 0$$

Void Space


$$\rho(\mathbf{u} \cdot \nabla \mathbf{u}) = -\nabla p + \mu \nabla^2 \mathbf{u}$$

$$\nabla \cdot \mathbf{u} = 0$$


Fluid velocity and pressure profiles

CFD model suggested different flow paths in the system

Flow Streaklines




Pressure Distribution / 1.0 ml/min



➤ Heterogenous flow distribution, asymmetrical

Multiple Flow Paths





Flow fraction:

Flow Path	% of Total Flow
Path 1	10.5
Path 2	36.0
Path 3	2.3
Path 4	51.2

Dye Binding Studies


- Cibacron Blue dye filtered through capsule
- Capsule cut-open to examine filter layers



Extend CFD model to predict protein breakthrough curves

Effect of module geometry on protein purification


Supracap™ 50 module



Stainless-steel module




(D)




Extend CFD model to predict protein breakthrough curves

Effect of module geometry on protein purification


Supracap™ 50 module




Stainless-steel module



Evaluate flow structure in both geometry




None uniform flow distribution



Uniform flow distribution

Protein Breakthrough

Dynamic binding capacity strongly influenced by flow distribution



- Langmuir binding kinetics used to account for protein binding

$$\frac{C_{eq}}{q_{eq}} = \left(\frac{1}{q_m} \right) C_{eq} + \left(\frac{1}{q_m} \right) \left(\frac{k_d}{k_a} \right)$$

- Stainless-steel module → sharper breakthrough curve with longer lag-time for initial breakthrough

Depth Filter Scale-Up

Pall Stax™




Area: 0.25 m²

Developing CFD model for large-scale depth filter


CFD model developed for flow and pressure profiles within the large-scale capsule

Stax module configuration





Inverted flow order

Streaklines in Stax module



Scale-up Effect: Flow Paths


Large-scale has much more uniform flow distribution




Flow Fractions		
	Stax	SupraCap
Path 1	0.4	10.5
Path 2	97.2	36.0
Path 3	0.01	2.3
Path 4	2.4	51.2

- Path 2 is desired flow path in both modules
- Stax module shows very little bypass (<2%) due to the small relative surface area for edge of Paths 3 and 4

Scale-Up Effects: Dye Binding Studies




Supracap (small-scale)



Stax module (large-scale)

Scale-Up Effects: Residence Time Distribution



- Residence time distribution evaluated using non-interacting tracer
- Supracap™ 50 module shows faster initial breakthrough due to bypass flow
- Stax™ module shows much sharper breakthrough curve due to more uniform flow distribution

Cell culture fluid (CCF) clarification

Small-Scale – PendoTech




Large-Scale – PendoTech




- CHO cell culture: 13×10^6 cells / mL, 60% viability
- Constant flux filtration: 0.6 and 1.2 L/m²/min

TMP profiles



- Pressure increases faster in the small-scale Supracap module due to non-uniformities in flow distribution
- Lower capacity at the higher filtrate flux → likely due to compressibility of deposited cells

DNA removal / breakthrough




- More rapid DNA breakthrough in Supracap module due to flow bypass
- Results validated that the flow path is more uniform in the large-scale module

Conclusions

- CFD model developed to describe flow profiles in both scale-down and pilot-scale depth filtration modules:
 - ✓ Significant bypass flow in small-scale module
 - ✓ Predictions validated using dye-binding and RTD studies
- Protein / DNA breakthrough curves evaluated by combining CFD analysis with Langmuir binding isotherm
 - ✓ Model predicts diffuse breakthrough in small-scale module
- Implications for scale-up
 - ✓ More uniform flow in pilot-scale module leads to greater capacity (based on TMP) and sharper breakthrough curves
 - ✓ General approach can lead to more accurate scale-up with reduced filter requirements

Acknowledgments


Prof. Andrew Zydney



Prof. Ali Borhan



Dr. Zydney's group



Merck colleagues



- Funding and assistance with the in-site experiments

