


cellicon® cell retention solution

Michale Brown
January 30th, 2023 AstraZeneca

Millipore
Sigma



MilliporeSigma is the U.S. and Canada Life
Science business of Merck KGaA, Darmstadt,
Germany.

Millipore®

Preparation, Separation,
Filtration & Monitoring Products

01 introducing the cellicon® cell retention solution!

Millipore®

Preparation, Separation,
Filtration & Monitoring Products

Cellicon® Cell Retention Solution

Cellicon® Filter Assembly

- Single-use, closed, gamma-irradiated assembly
- Flat-sheet tangential flow design
- 5 µm Durapore® membrane

Cellicon® Controller (lab scale)/ Mobius® Cell Retention System (process scale)

- Lab Scale/Process scale
- Allows for real-time process data
- Recirculation flow control, perfusate control (process scale)
- Intuitive software
- Remote connection capability

Initial Target Applications

- Upstream process intensification
- Monoclonal antibodies and r-proteins



Lab Scale

Process Scale

Scalable offering for perfusion up to 2000 L

Millipore
SIGMA

A complete solution
available across the
single-use bioreactor
size range for
perfusion
applications



CRS 50 L	CRS 200 L	CRS 500 L	CRS 1000 L	CRS 2000 L	
450 mm x 964 mm x 1588 mm	450 mm x 964 mm x 1588 mm	450 mm x 964 mm x 1588 mm	600 mm x 1064 mm x 1588 mm	700 mm x 1064 mm x 1588 mm	Dimensions (W x D x H)
(17.7" x 38.0" x 62.5")	(17.7" x 38.0" x 62.5")	(17.7" x 38.0" x 62.5")	(23.6" x 41.9" x 62.5")	(27.6" x 41.9" x 62.5")	

Bioreactor Size	3 L	50 L	200 L	500 L	1000 L	2000 L
Membrane Area	0.01m ²	0.2m ²	0.8m ²	1.9m ²	3.8m ²	7.6m ²
Levitronix® Pump	i30 SU	i100 SU	600 SU	600 SU	2000 SU	2 x 2000 SU
Connection Feed & Retentate Lines	SW	SW/AQ	SW/AQ	SW/AQ	SW/AQ	SW/AQ
Connection Perfusate Line	SW	SW	SW	SW	SW	SW
Recommended Crossflow Rate (L/min)	0.08	1.57	5.7	14	28	2 x 28
Recommended Max Flux (LMH)	20	20	20	20	20	20
Filter Hold-Up Volume (L)	0.014	0.24	0.73	1.7	3.5	7
Scale	PD/lab	Production				


SW: sterile welding
AQ: AseptiQuik®
sterile connector

Millipore®

Scalable offering for perfusion up to 2000 L

Millipore
SIGMA

A complete solution
available across the
single-use bioreactor
size range for
perfusion
applications



Bioreactor Size	3 L	50 L	200 L	500 L	1000 L	2000 L
Membrane Area	0.01m ²	0.2m ²	0.8m ²	1.9m ²	3.8m ²	7.6m ²
Levitronix® Pump	i30 SU	i100 SU	600 SU	600 SU	2000 SU	2 x 2000 SU
Connection Feed & Retentate Lines	SW	SW/AQ	SW/AQ	SW/AQ	SW/AQ	SW/AQ
Connection Perfusate Line	SW	SW	SW	SW	SW	SW
Recommended Crossflow Rate (L/min)	0.08	1.57	5.7	14	28	2 x 28
Recommended Max Flux (LMH)	20	20	20	20	20	20
Filter Hold-Up Volume (L)	0.014	0.24	0.73	1.7	3.5	7
Scale	PD/lab		Production			

SW: sterile welding
AQ: AseptiQuik®
sterile connector

Cellicon® Filter Assembly Family

Cellicon® Filter Assembly Contains:

- Flat sheet tangential flow filter designed for robust performance, easily scales across different process volumes
 - Lab scale: 1-3 L
 - Process scale: 50, 200, 500, 1000, and 2000 L (2 x 1000 L filters)
- Single-use components including:
 - Low shear levitating centrifugal feed pump
 - Pressure sensors for precise monitoring of filter performance
 - A perfusate flow sensor on process-scale assemblies for CSPR control
- Sterile connection methods to bioreactor:
 - Sterile weld (all sizes)
 - AseptiQuik® G connectors on feed and retentate (process-scale only)



Lab Scale



Process Scale



Cellicon® Filter Assembly Components

Process Scale

Cellicon® Filter

PendoTECH® SU
Pressure Sensors

C-Flex® 374
Weldable Tubing


AseptiQuik® G
Connectors

Bioprene® Tubing
for Perfusate
Pump

Levitronix® SU
Flowmeter


C-Flex® 374
Weldable Tubing

Levitronix® SU
Pump Head



Cellicon® Controller

Lab Scale





Specifications

- Control system designed to operate the lab-scale cell retention device
 - Single-use pump head and pressure sensors
 - Reusable clamp-on flow meter and pump motor
 - HMI control console with our custom user interface
- Scalability to commercial scales (3 – 2000 L bioreactors)
- Can be integrated into DeltaV™ Distributed Control System (DCS) and other DCS (via Modbus interface)


Mobius® Cell Retention Systems

Process Scale

- Control system designed to operate the process-scale cell retention device, containing:
 - a Levitronix® centrifugal feed pump drive
 - a non-invasive multi-use clamp-on ultrasonic feed flow meter
 - a Watson-Marlow peristaltic pump for perfusate
 - a connector for single-use ultrasonic flow sensor for perfusate
 - connectors for single-use pressure sensors in feed, retentate, and perfusate lines
- Supports external system connectivity over EtherNet/IP interface, allowing the data exchange with other systems
- Provided with our Bio4C ACE™ Software



Mobius® Cell Retention Systems Components



02 Application data


Millipore®

Preparation, Separation,
Filtration & Monitoring Products

Cellicon® Cell Retention Solution

A novel perfusion technology for process intensification


Millipore
SIGMA



Cellicon® Cell Retention Solution

A novel perfusion technology for process intensification



Millipore
SIGMA



High capacity

Cellicon® solution enables high viable cell densities in days

An N-1 perfusion run using CHO-S cell line was performed in a bioreactor across various scales with the Cellicon® Cell Retention Solution



Achieves high cell densities




Filter sized to meet the challenges of high VCD

High capacity

Cellicon® filter outperforms competitor device

N-1 Perfusion Runs:

	Membrane Area	Crossflow	Max Cell Density	Run Volume	Max Flux
Cellicon® filter (50 L)	1400 cm ²	0.925 LPM	51.5 e ⁶	50 L	24.2 LMH
Competitor	1300 cm ²	1.0 LPM	51.8 e ⁶	17 L	11.2 LMH





Cellicon® Cell Retention Solution achieves higher capacity and product yield

High capacity

Evaluation of Cellicon® Solution in N perfusion process

A steady-state perfusion run at 60×10^6 cells/mL using CHO-S cell line was evaluated in a bioreactor at 2 L working volume with the 3 L Cellicon® Cell Retention Solution.




- Completed 17-day process, 12 days of perfusion
- Achieved 2953 L/m^2 throughput at 2 psi TMP

Cellicon® Cell Retention Solution can support N perfusion applications


Evaluation of Cellicon® Solution in N process

A steady-state perfusion run at 100 e⁶ cells/mL using CHO-S cell line was evaluated in a bioreactor at 2 L working volume with the 3 L Cellicon® Cell Retention Solution.




Performance Data



Throughput vs. TMP



Protein Sieving




Performance Data



Cellicon® Cell Retention Solution

A novel perfusion technology for process intensification

Millipore
SIGMA



Scalable



Cellicon® Cell Retention Solution scales linearly from 3 to 2000 L



	PD	50	200	500	1000	2000
Max Working Volume (L)	2.7	50	200	500	1000	2000
Membrane Area (m ²)	0.01	0.2	0.8	1.9	3.8	7.6
Max Permeate Flow Rate (mL/min)	3.7	78.7	286.5	702	1404	2808
Crossflow Rate (LPM)	0.08	1.57	5.72	14	28	56
Levitronix® Pump Size	i30	i100	i600	i600	i2000	2xi2000
Wall shear	2467	2476	2484	2482	2482	2482

Key Scaling Parameters:


- **Channel height and length** are constant across all sizes
- **Membrane area** scaled in relation to bioreactor volume (37-43 cm²/L)
- **Crossflow rate** scaled by shear rate on membrane surface (2500 s⁻¹)
- Levitronix® **pump** selected based on RPM range

Scalable

Proven scalability using the Cellicon® Cell Retention Solution


A N-1 perfusion run using CHO-S cell line was evaluated in bioreactors at 50 & 200 L working volumes with the Cellicon® Cell Retention Solution

N-1 Process Cell Growth with Cellicon® Solution



Consistent cell growth and viability

N-1 Process Cellicon® Filter Performance (TMP vs. VCD)




No fouling and comparable performance

Cellicon® Cell Retention Solution

A novel perfusion technology for process intensification

Millipore
SIGMA




Robust

Filter provides additional capacity beyond N-1 process targets


A N-1 perfusion run using CHO-S cell line was evaluated in bioreactors at 50 & 200 L working volumes with the Cellicon® Cell Retention Solution

N-1 Process Cell Growth with Cellicon® Solution



High viability and cell growth across all scales

N-1 Process Cellicon® Filter Performance (TMP vs. VCD)




Filter appropriately sized for application

Robust


Consistent performance with different cell lines

Multiple N-1 perfusion runs were performed using CHO-S and CHOZN® & UCOE® Combined Platform cell lines evaluated in bioreactors at 2-2.5 L working volume targeting 20pL/cell/day CSPR with the 3 L Cellicon® filter to evaluate performance.

Viable Cell Density and Viability



Throughput vs TMP




High viability and cell growth across all cell lines

Cellicon® Cell Retention Solution

A novel perfusion technology for process intensification

Millipore
SIGMA



Low shear

Crossflow rates determined to match shear rate across scales



Cellicon® filter assembly	3 L	50 L	200 L	500 L	1000 L	2000 L
Flow Rate (L/min)	0.08	1.57	5.72	14.00	28.00	56.00
Number of Channels	1	3	11	27	54	108
Flow Velocity (cm/s)	12.5	14.7	14.7	14.7	14.7	14.7
Wall Shear (1/s)	2467	2476	2484	2482	2482	2482
Filter Area (m ²)	0.01	0.2	0.8	1.9	3.8	7.6

$$\text{Wall shear} = \frac{6v}{h}$$


v = flow velocity
= Q_L /flow area
 h = channel height

Shear maintained across scales ensuring consistent performance


Low shear

Shear rate effect on cell culture performance

Cell Density, Viability, and Doubling Time



LDH Profile




Low shear in the flow path preserves cell health

Cellicon® Cell Retention Solution

A novel perfusion technology for process intensification

Millipore
SIGMA



Complete monitoring and control


Bio4C ACE™ Software



Navigation Bar

Banner

Process Monitor and Control




Run Panel

Simplified Trends


Complete monitoring and control

Integrated instruments ensures complete process control


50 L Pressure Data



50 L Feed Flow & Speed



50 L Perfusate Flow & Speed




A N-1 perfusion run using CHO-S cell line was evaluated in a bioreactor at 50 L working volume with the Cellicon® Cell Retention Solution

Complete monitoring and control delivers robust performance

Cellicon® Cell Retention Solution

A novel perfusion technology for process intensification

Millipore
SIGMA



Simple & Quick

Ready to set-up and process in minutes

1

System installation:

- Minimal utility requirements, only needs power supply

2

Filter installation:


- Ships as a complete assembly (filter, pressure sensors, pump head, and flow meter)
- Gamma irradiated
- Sterile welding and sterile connector
- No need to pre-wet/pre-flush the membrane

Ready to use in minutes

Cellicon® Cell Retention Solution

A novel perfusion technology for process intensification

Millipore
SIGMA



03 case study

Millipore®


Preparation, Separation,
Filtration & Monitoring Products

Case study


Experiment designed to evaluate an intensified upstream process

Experimental Condition: Perfused N-1 Seed Train

Transfer at ~ 50 million cells/mL





Control Condition: Conventional Seed Train



Case study

Achieves high cell density without challenging filter


Bioreactor: 2L working volume
Media: Cellvento® 4CHO-X COMP expansion media, target 20 pL/cell/day CSPR
Cell Retention: Bench-scale Cellicon® Cell Retention Solution at 80 mL/min crossflow
Inoculation Density: 0.2e⁶ cells/mL
CHO-S cells transferred on day 8 to fed-batch production bioreactor (~55 e⁶ cells/mL)




Case study

N bioreactors were seeded at conventional and high cell densities

Experimental Condition: High Seed Intensified Fed-Batch





Historical Control Condition: Conventional Fed-Batch



Case study



High seed condition achieves peak VCD faster



- ✓ Peak VCD achieved 4 days earlier for high seed production bioreactor with higher peak VCD observed
- ✓ Similar trends in viability across conditions and compared to historical fed-batch runs

Case study



Final protein concentrations achieved four days sooner with high seed perfused N-1 condition



- ✓ Similar final titer between intensified process and control or standard seed perfused N-1 was achieved four days earlier in high seed condition
- ✓ Comparable specific productivity in high seed condition, low seed, and historical data
- ✓ Utilizing Cellicon® Cell Retention Solution did not impact productivity

Case study


Product quality is maintained during intensified processing



- ✓ SEC analysis shows consistent product quality across conditions
- ✓ Consistent charge profiles observed between conditions
 - No impact between high and low seed densities
 - Comparison to historical control is within the normal range of variability for this process

Case study

Product quality is maintained during intensified processing



- ✓ Process intensification strategy did not impact the product glycan profile when compared to the control process

Case study

Intensified processing benefits

Seed Train



Less consumables and equipment needed to complete the seed train

- Reduced flasks and/or bioreactors needed
- Reduced risk of contamination with fewer process steps



Increased flexibility within seed train

- Enables high seed inoculations, multiple or larger bioreactors
 - Achieved with higher cell concentrations at lower volumes
- Reduced time or increased productivity in production process

Production



Save time and labor

- Reduced production time with high seed conditions
- Allows for more batches per year



Intensified processing benefits **amplify at large-scale production**

- Cost and time savings
- Reduced facility requirements – less space needed or increased manufacturing capabilities

Thank You!

